

Instructions for Use

HYCOAT®

Topical Wound Management Solution

Hyaluronate sodium sterile solution

Sodium hyaluronate/hyaluronic acid (HA) plays an important role in the wound healing process. Sodium hyaluronate acts as a lubricant on the skin by restoring moisture and aiding in prevention of friction or abrasion. It is a polysaccharide that is distributed widely in the extracellular matrix of connective tissue. HYCOAT® helps manage surgically induced and chronic wounds by providing a moisture rich environment and is easy to handle and deliver.

Description:

HYCOAT® is available in a 10 mL solution containing 50 mg, and a 2mL solution containing 20mg of sodium hyaluronate having a pH between 6.5-7.25 and is preservative free.

Indications:

HYCOAT® is indicated for use in the topical management of chronic and acute wounds including:

- ❖ Pressure ulcers (Stages I-IV)
- ❖ Acute wounds (such as cuts, abrasions, and lacerations)
- ❖ First and Second Degree Burns
- ❖ Surgical wounds (including donor sites)
- ❖ Venous stasis ulcers
- ❖ Autograft procedures

Precautions and Contraindications:

Discard unused portion. HYCOAT® is to be dispensed by prescription only.

No known side effects. Not intended for use as a long-term or permanent dressing.

Product Administration:

Using HYCOAT®:

1. Cleanse wound/surgical site.
2. Withdraw the HA solution.
3. Dispense topically onto the wound site.
4. Cover wound area with a non-adherent dressing.
5. Change dressing as needed, repeating application of HYCOAT®. *Infected wounds are treated per physician direction.

Storage:

- ❖ Store at temperatures between (15° – 30° C); Protect from freezing; Sterile when unopened.

HYCOAT® is manufactured in the USA for Hymed® and is a registered trademark of The Hymed Group Corporation.

The information given herein is believed to be reliable. However, no guarantee is made or liability assumed. Hymed® does not know all of the uses to which its products may be utilized or the condition of use. Therefore Hymed® makes no claims or warranty concerning the fitness or suitability of the product for a particular use or purpose other than those outlined.